

# INTELLIGENT FACE MASSAGER

## User Manual & Product Specifications

### PRODUCT PARAMETERS

|                 |                     |
|-----------------|---------------------|
| Product Name    | Smart Face Massager |
| Model           | LM-681              |
| Adapter Voltage | 5V = 1A             |
| Rated Voltage   | 3.7V =              |
| Nominal Power   | 4W                  |

### CORE FACE-LIFT MODES

- Automatic Mode
- V-Face Mode
- Compact Mode
- Pull-Up Mode
- Shaping Mode
- Micro-current Pulse

### EXQUISITE V-FACE LITTLE SECRET (INTENSITY ADJUSTMENT)

For individuals with different facial needs, 12 levels of intensity adjustment are designed according to differences in facial skin and thickness quality.

**Levels 1–3:** Targeted for first-time users or individuals with relatively sensitive skin, thin stratum corneum, and fine skin texture.

**Levels 4–5:** Recommended for users who have used the product for 3–5 weeks, targeting people with moderate masseter muscle hypertrophy and a double chin.

**Levels 6–8:** Targeted for people who have used the product for about 3 months, or those with a wide and flat jawline under severe stress.

**Levels 9–12:** Targeted for advanced use (half a year or more), reinforcing and consolidating the face-shining effects.

## INSTRUCTIONS FOR USE

- 1 Before using the product, thoroughly clean your face and ensure a good application of skincare and moisturizing products.
- 2 Wet the massage sheet before use. **(It is strictly forbidden for liquid skincare products such as water, milk, or cream to infiltrate into the interior of the device to avoid affecting the product's lifespan!)**
- 3 Place the massage sheet against the masseter muscle, adjust the position, pull the bandage upwards, and secure it. The bandage can be adjusted according to the size of the head.
- 4 **Button Controls:**
  - ① Power Button
  - ② Increase Intensity Key
  - ③ Decrease Intensity Key
  - ④ Mode Cycle Key

### IMPORTANT TIPS & MAINTENANCE

1. **Skincare Application:** Please clean your face before use and properly apply moisturizing skincare products. **It is strictly forbidden to let liquid skincare products like water, milk, or cream penetrate the interior of the device shell to prevent short circuits and ensure a long lifespan!**
2. **Sheet Preparation:** Dampen the massage sheet slightly before application to adjust the elasticity of the bandage, fit the product fully to the face, and adjust the massage intensity gradually from low to high.
3. **Charging Guidelines:** Turn off the product when charging. **It is strictly prohibited to use the device while it is charging!** The red charging symbol flashes during the charging process. When charging is complete, the charging symbol turns solid green. After a full charge, you must pull the plug out in time to avoid prolonged charging and shortening the product's battery life.
4. **Auto Shut-off:** The product will automatically shut down after 15 minutes of continuous operation.