

3D LIGHT THERAPY ANTI-AGING MASK

USER MANUAL • MODEL: V-F2509

Before use, please read this manual in its entirety. Welcome to using this product.

This manual provides detailed information on the product's features, how to use it, and precautions. Please carefully read this manual before using the device. After reading, please keep it properly for future reference. The images in this manual are for reference only, please refer to the actual product.

Product Introduction

This product is a household electronic beauty device that uses green light, yellow light, red light, and near-infrared light. It applies the effects of different wavelengths in the four spectra to provide beauty and repair for various skin types. The 3D light therapy anti-aging mask is a non-invasive and painless skin care method without chemical additives or light pollution, and it is highly portable and effective.

Product Specifications

Product Name	3D Light Therapy Anti-Aging Mask	LED Lamp Beads	328 PCS
Model	V-F2509	Timing	5 / 10 / 15 MIN
Product Size	L220 × W169 × H119 mm	Input	5V ... 1A
Net Weight	Mask: 290 ± 10g Controller: 90g	Charging Time	3.5 Hours
Battery Capacity	2000mAh, 3.7V ...	Operation Temp.	0 – 40°C
Skin Type	All Skin Types	Wavelengths	Green: 525nm Yellow: 590nm Red: 660nm Near-Infrared: 850nm
Package Includes	Mask ×1, Charging Cable ×1, Controller ×1, User Manual ×1, Male-female Velcro ×1 (This accessory can be used when the bonding length is insufficient)		

Product Interface & Controls

Controller Features

- **Time Display:** Shows remaining session time.
- **Charging Indicator:** Displays the power and charging status.
- **Mode Button (M):** Short press to switch between the 4 care modes.
- **Timer Button (T):** Short press to adjust the session time (5, 10, or 15 minutes).
- **Power Button (⏻):** Long press for 1.5 seconds to turn the device on or off.

Note: The Type-C port on the controller is compatible with both input charging and output connection to the mask. The internal smart chip automatically detects the required mode.

Care Steps

Step 1: Makeup Removal and Cleansing

Before using the device, it is essential to thoroughly remove makeup and cleanse the skin to ensure a clean face for light treatment.

Step 2: Connecting the Mask

Connect the 3D light therapy anti-aging mask to the controller using the Type-C cable. Check to ensure that the device powers on properly, verify the battery level, and ensure all lights are functioning correctly. Adjust the Velcro straps on the mask to achieve a secure and comfortable fit on your face.

Step 3: Activate Care Mode

Press and hold the power button for 1.5 seconds to turn on the device. Use the **T button** to adjust the timer (5/10/15 minutes), and press the **M button** to switch light modes and start your beauty care routine.

Usage Optimization Tips:

1. Before first use, shine the device on your face for one minute. If there are no adverse reactions, you can proceed with regular use.
2. Start with the default 5-minute setting. Once your skin adapts, adjust the time according to your skin's tolerance and specific needs.
3. Avoid excessive light exposure time and frequency. Always adjust the duration based on your skin's actual condition.
4. It is recommended to use the product 3–4 times a week for 10–15 minutes per session.
5. Temporary faint imprints around the eyes due to tight wear are completely normal and will fade.

Step 4: Hydrate After Treatment

After the photon treatment, you can use a hydrating face mask for extra moisture, or simply follow your regular skincare routine for enhanced benefits. Prioritize deep hydration and sun protection in your daily skincare routine.

Care Modes Detail

Mode 1: Calm, Soothe and Cleanse Acne

Wavelength: 525nm Green Light

Benefits: Green light provides anti-inflammatory, antibacterial, soothing, and calming effects. It cleanses the skin and effectively inhibits excessive sebum production.

Mode 2: Brighten

Wavelength: 590nm Yellow Light

Benefits: Yellow light enhances cellular oxygen alternation, whitens, lightens pigmentation, and improves microcirculation. It refines rough skin textures and addresses poor circulation.

Mode 3: Rejuvenation

Wavelength: 660nm Red Light

Benefits: Red light enhances blood circulation, accelerates metabolism, and stimulates collagen production. This effectively reduces fine lines and reveals radiant, glowing, and rejuvenated skin.

Mode 4: Tighten & Firming Skin

Wavelength: 660nm Red Light + 850nm Near-Infrared Light

Benefits: Red and near-infrared light penetrates deeply into the tissue, stimulating cell membranes and accelerating cellular self-repair to firm fine lines and loose skin.

Charging Method

- This product is equipped with a versatile bidirectional Type-C port, integrating both charging and power supply functions.
- Before first use, it is highly recommended to fully charge the controller.
- For charging, please use a 100-240V AC universal **5V ... 1A** adapter (Note: Adapter is not included).
- The charging indicator light flashes during the charging process and remains solid when charging is fully complete.

Care, Storage & Precautions

Maintenance and Storage

- After care, place the device back into its packaging box.
- If dirty, wipe the surface clean with a soft, damp cloth. Do not rinse or submerge the device in water.
- Keep in a cool, dry place away from moisture and direct sunlight. Store out of reach of children.

Important Precautions

- **Not suitable for children.**

- Pregnant or breastfeeding women should use this device with caution.
- People with epilepsy should use with caution, as certain light waves and flashing frequencies may trigger seizures.
- Do not use if you have known light sensitivity (photosensitivity).
- Individuals undergoing professional dermatological treatment or in the acute phase of a skin condition (such as lupus erythematosus, photosensitive eczema, vitiligo, thyroid disease, etc.) are advised not to use this product except under strict medical guidance.
- If you feel any discomfort or unusual irritation during exposure, stop the operation immediately and consult a medical professional.

Warranty & After-Sales

Thank you for choosing our products. We are pleased to offer you a **365-day free warranty service**. If you have any inquiries, issues, or concerns, please don't hesitate to reach out to our customer support team.